

Laboratory Chest Freezer ProfiLine Taurus PLTA


- Static cooling and manual defrost Most reliable and less complicated compressor technology, low-noise
- Micro-processor driven temperature controller with digital display
- Visual/audible alarm device
 - Alarm warning for excess and insufficient temperature
 - No-volt remote alarm contact
- External housing made of white, galvanized steel plate
- Interior made of white plastic
- Lockable lid and spring-loaded hinges to allow opening at any angle
- With castors, supplied separately (except for PLTA 0962)
- Storage baskets optional (see below)



Model ProfiLine Taurus	PLTA 0962	PLTA 1462	PLTA 3362	PLTA 4362	PLTA 5562
Temperature range	°C:-30/-60	-30/-60	-30/-60	-30/-60	-30/-60
Temperature range -40 °C down to -86 °C see separate data sheet					
Content litres (gross):	74	133	296	383	532
Ext. dimensions (WxDxH) mm:	555x635x885	720x695x890	1260x695x890	1560x695x890	1705x655x865
Int. dimensions 1 (WxDxH) mm:	387x387x470	350x435x650	880x435x650	1190x435x650	1340x490x700
Int. dimensions 2 (WxDxH) mm:	—	200x435x400	200x435x400	200x435x400	200x490x450
Height at open top mm:	1480	1500	1500	1500	1540
Baskets optional (max. quantity)	: 4 (Typ A)	4 (Typ B)	13 (Typ B)	16 (Typ B)	19 (Typ C)
Weight kg:	53	55	76	88	107
Mains connection	: < 1/N/PE ~ 230 V 50 Hz >				
Shipping dimensions (WxDxH) mm:	680x640x1160	810x760x1090	1360x760x1060	1660x760x1060	1720x800x1150
Shipping weight kg:	83	85	114	132	131

Options / Accessories

Options / Accessories	Cat. No.
Power failure alarm	See models PLTA ..63 on page 5.2
ProfiTec Temperature Monitoring and Recording System (see separate brochure). USB data logger up to -200 °C with LCD display.	GPTEC4A ZXX2-045-010
ProfiCaller telephone dialler.	On request ZXX2-020-2L
Pt-100 temperature sensor, additionally installed for external temperature recording.	ZXX2-093
CO2 back-up system, free-standing.	ZPLTA2-094
Facilities for customer connection of CO2 back-up system.	ZPLTA5-21-09xx
Transport carriage with 4 swivel castors (for PLTA 0962).	on request
4 castors, 2 lockable (for PLTA 5562).	on request
Storage basket (type A), 192x384x190 mm (WxDxH).	ZXX4-084
Storage basket (type B), 220x430x200 mm (WxDxH).	ZXX4-085
Storage basket (type C), 220x480x215 mm (WxDxH).	ZXX4-093




ProfiTec Intelligent Temperature Monitoring and Recording System


- For up to 4 refrigerators and/or freezers
- Readings are stored over a period up to 781 days
- Warning for excess and insufficient temperature
- Includes No-volt remote alarm contact
- Can be fed to PC via RS 232 port



Storage sets with matching vertical racks made of stainless steel




Number of racks when fully equipped



Type of box
ALPHA Box
BRAVO Box
EPPi® Box



Type of box
ALPHA Box
BRAVO Box
EPPi® Box



Type of box
ALPHA Box
BRAVO Box
EPPi® Box



4 or 5

Type of box	Cat. No.
ALPHA Box	TA50-08AN
BRAVO Box	TE50-08AN
EPPi® Box	TE50-08AN

Type of box	Cat. No.
ALPHA Box	TA75-05AN
BRAVO Box	TB75-05AN
EPPi® Box	TE80-05AN

Type of box	Cat. No.
ALPHA Box	TA100-04AN
BRAVO Box	TE102-04AN
EPPi® Box	TE102-04AN

PLTA 1462 on request



18

Type of box	Cat. No.
ALPHA Box	TA50-11AN
BRAVO Box	TE50-11AN
EPPi® Box	TE50-11AN

Type of box	Cat. No.
ALPHA Box	TA75-07AN
BRAVO Box	TB75-07AN
EPPi® Box	TE80-07AN


Type of box	Cat. No.
ALPHA Box	TA100-03AN
BRAVO Box	TE102-03AN
EPPi® Box	TE102-03AN

3

Type of box	Cat. No.
ALPHA Box	TA50-06AN
BRAVO Box	TE50-06AN
EPPi® Box	TE50-06AN

Type of box	Cat. No.
ALPHA Box	TA75-04AN
BRAVO Box	TB75-04AN
EPPi® Box	TE80-04AN

Type of box	Cat. No.
ALPHA Box	TA100-03AN
BRAVO Box	TE102-03AN
EPPi® Box	TE102-03AN



24

Type of box	Cat. No.
ALPHA Box	TA50-11AN
BRAVO Box	TE50-11AN
EPPi® Box	TE50-11AN

Type of box	Cat. No.
ALPHA Box	TA75-07AN
BRAVO Box	TB75-07AN
EPPi® Box	TE80-07AN

Type of box	Cat. No.
ALPHA Box	TA100-03AN
BRAVO Box	TE102-03AN
EPPi® Box	TE102-03AN

3

Type of box	Cat. No.
ALPHA Box	TA50-06AN
BRAVO Box	TE50-06AN
EPPi® Box	TE50-06AN

Type of box	Cat. No.
ALPHA Box	TA75-04AN
BRAVO Box	TB75-04AN
EPPi® Box	TE80-04AN

Type of box	Cat. No.
ALPHA Box	TA100-03AN
BRAVO Box	TE102-03AN
EPPi® Box	TE102-03AN

PLTA 5562 on request