



Model No.	ML0801WU	ML1501WU	ML3001WU
Temperature range	+1°C to +5°C	+1°C to +5°C	+1°C to +5°C
Forced air circulation	Yes	Yes	Yes
Defrost (automatic / manual)	automatic	automatic	automatic
Electronic controller	Yes	Yes	Yes
Digital display	Yes	Yes	Yes
Temperature alarm visual / audible	Yes	Yes	Yes
Excess and insufficient temp. alarm	Yes	Yes	Yes
Open door alarm	No	Yes	Yes
No-volt remote alarm	No	No	No
External finish	painted steel plate	painted steel plate	painted steel plate
Colour	white	white	white
Lock	Yes	Yes	Yes
Door type	solid door	solid door	solid door
Door hinges	changeable	changeable	changeable
Door seal	not changeable	not changeable	not changeable
Internal finish	plastic, white	plastic, white	plastic, white
Interior light	No	No	Yes
Spark free	No	No	No
Shelves (number)	2	3	5
Shelf positions (number)	2	6	12
Shelf size (mm) (WxDxH)	440x275x6	520x320x6	497x340x6
Drawers (number)	0	0	1
Wire baskets (number)	1	1	0
Door bins (number)	0	2 (fixed)	4 (adjustable)
Optional drawer system	Yes	Yes	Yes
Adjustable feet (number & position)	No	2 front feet 2 rear rollers	2 front feet 2 rear rollers
Gross capacity (l)	82	151	288
Net capacity (l)	67	145	284
External width (mm)	502	635	595
External depth (mm)	540	595	670
Door open depth (mm)	976	1165	1185
External height (mm)	660	845	1570
Internal dimensions (mm) (WxDxH)	419x315x540 step (DxH): 165x215	500x450x705 step (DxH): 165x210	480x445x1382 step (DxH): 200x260
Product weight (kg)	27.5	37	59
Voltage supply	1/N/PE-230V 50Hz	1/N/PE-230V 50Hz	1/N/PE-230V 50Hz
Packed width (mm)	700	700	700
Packed depth (mm)	750	750	750
Packed height (mm)	825	1022	1735
Packed weight (kg)	34.5	44	66
Refrigerant type CFC / HCFC free	R134a	R600a	R600a



ML0801WU



ML0801WU with Inventory System (optional)



ML1501WU



ML3001WU



ML3001WU with Inventory System (optional)



Recommendations for additional storage sets can be found in our separate catalogue. Please ask for a personal fitting option.

