



ML0806WN



ML0806WN with Inventory System (optional)



ML1006WN



ML1006WN with Inventory System (optional)



ML1506WN



ML1506WN with Inventory System (optional)



ML3006WN

Model No.	ML0806WN	ML1006WN	ML1506WN	ML3006WN
Temperature range	+2°C to +10°C	+2°C to +10°C	+2°C to +10°C	+2°C to +10°C
Forced air circulation	No	No	No	No
Defrost (automatic / manual)	automatic	automatic	automatic	automatic
Electronic controller	Yes	Yes	Yes	Yes
Digital display	Yes	Yes	Yes	Yes
Temperature alarm visual / audible	Yes	Yes	Yes	Yes
Excess and insufficient temp. alarm	Yes	Yes	Yes	Yes
Open door alarm	No	No	Yes	Yes
No-volt remote alarm	No	No	Yes	Yes
External finish	painted steel plate	painted steel plate	painted steel plate	painted steel plate
Colour	white	white	white	white
Lock	Yes	Yes	Yes	Yes
Door type	solid door	solid door	solid door	solid door
Door hinges	changeable	changeable	changeable	changeable
Door seal	not changeable	not changeable	not changeable	not changeable
Internal finish	plastic, white	plastic, white	plastic, white	plastic, white
Interior light	No	No	No	No
Spark free	Yes	Yes	Yes	Yes
Shelves (number)	2	2	3	5
Shelf positions (number)	2	4	6	12
Shelf size (mm) (WxDxH)	438x274x6	427x374x6	520x320x6	497x340x6
Drawers (number)	0	0	0	1
Wire baskets (number)	1	1	1	0
Adjustable door bins (number)	0	0	2	4
Optional drawer system	Yes	Yes	Yes	Yes
Adjustable feet (number & position)	No	2 front feet	2 front feet 2 rear rollers	2 front feet 2 rear rollers
Gross capacity (l)	82	107	151	288
Net capacity (l)	67	103	145	284
External width (mm)	502	495	595	595
External depth (mm)	540	620	595	634
Door open depth (mm)	976	1050	1165	1185
External height (mm)	660	830	845	1570
Internal dimensions (mm) (WxDxH)	415x340x540	403x400x540	500x450x705 step (DxH): 165x185	480x445x1382 step (DxH): 175x255
Product weight (kg)	27.5	31.5	36	58
Voltage supply	1/N/PE-230V 50Hz	1/N/PE-230V 50Hz	1/N/PE-230V 50Hz	1/N/PE-230V 50Hz
Packed width (mm)	700	700	700	700
Packed depth (mm)	750	750	750	750
Packed height (mm)	825	1000	1022	1735
Packed weight (kg)	34.5	38.5	43	65
Refrigerant type CFC / HCFC free	R134a	R134a	R600a	R600a



Recommendations for additional storage sets can be found in our separate catalogue. Please ask for a personal fitting option.

